<div><div><span><span> </span></span></div><div><b>BAAT Medical Products B.V.</b></div></div> <p>Fabrikantenadres:      7555 RJ Hengelo      The Netherlands   </p> <p>Phone: +31 6088 565 60      www.baatmedical.com</p>	<div><div><span><span> </span></span></div><div><b>BAAT medical</b></div><div><b>PRODUCTS</b></div></div>
--	---

**CE** **CE** 0344

FU-PR00110\_nvnl\_07\_2015

## ENGLISH

### Important Medical Information

PRE-OPERATIVE PRECAUTIONS

This package insert includes information on basic precautions and reprocessing instructions for reusable medical devices from BAAT Medical Products. Please consult the applicable surgical technique for selection of a device and check the full labelling for other important information. Surgical technique brochures may be requested from a distributor or from BAAT Medical Products directly. Those using brochures published more than two years before the surgical intervention are advised to obtain an updated version.

BAAT Medical Products devices can only be used by persons who are fully familiar with the surgical technique required and who have been trained in this role. The operating surgeon must take care not to use the instruments to exert inappropriate stress on the patient or the implants and must comply completely with the operating procedure described in the surgical technique brochure. For example, the forces exerted when repositioning an instrument in-situ must not be excessive as this is likely to cause injury to the patient.

To reduce the risk of brucella, care must be taken to not distend the implants or, risk to not score them with the instruments unless otherwise specified by the applicable surgical technique. Extra care must be taken when the devices are used near vital organs, nerves or vessels. Unless otherwise specified on the label, the devices may be reused after decontamination, cleaning and sterilization.

WARNINGS (U.S.A.)
   
Federal law restricts this device to sale by or on the order of a licensed physician

CONTACT:
   
See brochures for telephone and address of local representative/distributor, or contact BAAT Medical Products.

**(Re)processing instructions for reusable medical devices**
  
Devices which are delivered sterile/sterile are available in a minimum quantity of 25k/g of gamma radiation from a cobalt 60 source or have been sterilized by a dynamic-air removal steam sterilization method (see product label and sterilization method applied).

Devices not labelled as sterile or non-sterile. The packaging of all sterile devices should be inspected for flaws in the sterile barrier or expiration of shelf life before opening. In the presence of such a flaw or expiration of shelf life, the product must be assumed non-sterile.

In the event of contamination of reusable devices, or expiration of shelf life in the product use supplied non-sterile, the device must be subjected to an appropriate and validated cleaning and sterilization procedure. Use the instructions chart below for reprocessing medical devices.

In the event of contamination of a single use device with body fluids, the device is considered unclean, and must be discarded.

Reusable device	Supplier not sterile	Expiration shelf life	Contaminated with body fluids	Contaminated after use
Response	Response	Response	Response	Response
Single use device	Response	Response	Discard	Response

Table 1 - (Re)processing or discard of contaminated devices

**CAUTIONS**

- Long, narrow cannulae and blind holes require particular attention during decontamination.
- Always clean cannulated devices intra-operatively to prevent accumulation of debris in the cannulation.
- Mechanical instruments should be cleaned separately from other instruments.
- For good maintenance, reusable instruments should be pre-cleaned, cleaned and sterilized immediately after surgery.

**Limitations on reprocessing**

- Reprocessed packaging has minimal effect on these products. End of life is normally determined by wear and damage due to use.

<b>Preparation at the point of use</b>	<b>INSTRUCTIONS CHART FOR (RE)PROCESSING MEDICAL DEVICES</b>	
	Pre-cleaning	Pre-sterilization
<b>Contaminated and transportation</b>	<ul style="list-style-type: none"><li>Remove gross soil</li></ul>	<ul style="list-style-type: none"><li>It is recommended that the devices are reprocessed as soon as it is reasonably practical following use</li> <li>Disassemble when possible following assembly/Disassemble instructions</li></ul>
	Preparation for cleaning	Pre-sterilization (method preference)
<b>Cleaning – automated</b>	<ul style="list-style-type: none"><li>Pre-cleaning</li> <li>Reuse excess soil from device.</li> <li>Soak in cleaning solution for 10 minutes</li> <li>Use noncaustic brush</li> <li>Reuse with demineralized water</li></ul>	<ul style="list-style-type: none"><li>Apply a suitable quantity of surfactant gel solution to all surfaces ensuring that hinged instruments are cleaned in both open and closed positions</li> <li>Pay close attention to threads and hard to reach areas.</li> <li>Clean cannulae and holes using an appropriate brush ensuring that full depth of the feature is reached</li> <li>Reuse with demineralized water.</li></ul>
	Washer-disinfectant	<ul style="list-style-type: none"><li>Washer-disinfectant</li> <li>Load devices such that hinges are open and cannulae and holes can drain</li> <li>Rin cycle at minimum 92 °C for 10 minutes</li> <li>When unloading, check cannulation, holes, etc. for complete removal of visible soil.</li> <li>Use necessary repeat cycles or use manual cleaning</li></ul>
Washing	<ul style="list-style-type: none"><li>Pre-cleaning</li> <li>Reuse excess soil from device.</li> <li>Soak in cleaning solution for 10 minutes</li> <li>Use noncaustic brush</li> <li>Reuse with demineralized water</li></ul>	<ul style="list-style-type: none"><li>Pre-cleaning</li> <li>Reuse excess soil from device.</li> <li>Soak in cleaning solution for 10 minutes</li> <li>Use noncaustic brush</li> <li>Reuse with demineralized water</li></ul>

Manual cleaning

- Soak in ultrasonic bath for 15 minutes in mild detergent
- Use noncaustic brush, apply detergent solution to all surfaces ensuring that hinged instruments are cleaned in both open and closed positions
- Pay close attention to threads and hard to reach areas.
- Clean cannulae and holes using an appropriate brush ensuring that full depth of the feature is reached
- Reuse with demineralized water.

Drying

- When drying is achieved as part of a washer-disinfecter cycle, do not exceed 120 °C
- Dry devices until completely dry
- When drying is achieved as part of a washer-disinfecter cycle, do not exceed 120 °C

Maintenance

- Apply a small quantity of surgical grade lubricant oil to hinges and threaded sections.
- Inspect and, where necessary, clean damaged device.

**Inspection and testing**

- Hinged instruments: Check for smooth movement of hinge without excessive "play".
- Locking and ratchet mechanisms should be checked for function.
- Threaded devices must run smoothly.
- All devices: Visual inspection for damage and wear. Cutting edges should be free of nicks and present a continuous edge.

Check devices which are long, slender features (particularly rotating instruments) for distortion. When devices form part of a larger assembly, check assembly with mating components.
   
Assemble where possible following assembly/disassembly instructions.
   
Single A standard packaging material may be used. Ensure that the pack is large enough to contain the device without crushing the metal.
   
If sets, devices may be loaded in dedicated trays or general purpose sterilization trays. Ensure that cutting edges are protected. Double wrap the trays using appropriate method.

**Sterilization (Preferred method)**

- Dynamic-air removal steam sterilization method (pre-vacuum)
- Exposure time: 3 minutes
- Temperature: 134 °C (274 °F)
- Pressure: 3.0 bar (45 psi)
- Drying time: 16 minutes
- Control moisture
- Control storage time

**Drying**

- When drying is achieved as part of a washer-disinfecter cycle, do not exceed 120 °C
- Dry devices until completely dry
- When drying is achieved as part of a washer-disinfecter cycle, do not exceed 120 °C

**Maintenance**

- Apply a small quantity of surgical grade lubricant oil to hinges and threaded sections.
- Inspect and, where necessary, clean damaged device.

Disinfection method	Lidless not sterile	Overridden GHS/Inform	Contaminated with liquid contents	Contaminated with other contents
Method/s preferred	Response	Response	Response	Response
Product/Disinfectant	Response	Response	Response	Response

Table 1 - Disinfection of overriden or disinfected contaminated medical devices

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Table 1 - (Re)processing or disinfection of disipositi contaminated

PRECAUTIONS	Long narrow cannulae and blind holes require particular attention during decontamination.
Reinigungs – manuell	Pre-cleaning
Dispositiv reutilizable	Formio not sterile
Dispositiv reutilizable	Response
Dispositiv reutilizable	Response

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**Entretien**

- Apply a suitable quantity of surfactant gel solution to all surfaces ensuring that hinged instruments are cleaned in both open and closed positions.
- Pay close attention to threads and hard to reach areas.
- Clean cannulae and holes using an appropriate brush ensuring that full depth of the feature is reached.
- Reuse with demineralized water.

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**AMM Technical Information Report (TR) 12**
  
EN 285 Sterilization - Steam sterilizers - Large sterilizers
   
EN 15883 Washer-Disinfectors
   
EN 868 Packaging for terminally sterilized medical devices
   
ISO 11607 Packaging for terminally sterilized medical devices

**EXPLANATION OF SYMBOLS**
  
EN 980 Symbols for use in the labelling of medical devices

Dispositivo reutilizable	Supplier not sterile	Overridden GHS/Inform	Contaminated with liquid contents	Contaminated with other contents
Response	Response	Response	Response	Response
Dispositivo reutilizable	Response	Response	Response	Response

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